


SUPPLEMENTARY MATERIAL

Figure S1. Basophil activation test flow cytometry plots



A) BAT plots in the patient with the positive control N-formylmethionine-leucyl-phenylalanine (fMLP). Basophil were selected as CD203c+ cells. This test is considered positive if basophil degranulation is >5%. B) BAT plots for rocuronium, sugammadex and rocuronium-sugammadex mixture (SR-M) complex. This data shows no basophil response for rocuronium nor sugammadex exposure, but reflects a positive result to the SR-M complex.