

SUPPLEMENTARY MATERIAL




Figure 1S: A) Flow cytometric gating strategy used for the basophil activation test (BAT). Basophils were gated as CD123⁺/HLA-DR⁻ cells (first and second columns), and CD63 was used as a marker of basophil degranulation. B) CD63 basophil activation in response to fMLP, Fc ϵ R1 (positive controls) and 50 $\mu\text{g/ml}$ paclitaxel stimulation in a paclitaxel-allergic patient. C) CD63 basophil activation in response to fMLP, Fc ϵ R1 (positive controls) and 50 $\mu\text{g/ml}$ paclitaxel stimulation in healthy controls.